

## Product datasheet for **HK203730**

### Signal Peptide Peptidase (HM13) Human qPCR Template Standard (NM\_030789)

#### Product data:

|                  |   |
|------------------|---|
| Product Type:    | qPCR Template Standards   |
| ACCN:            | NM_030789   |
| Synonyms:        | H13; IMP1; IMPAS; IMPAS-1; MSTP086; PSENL3; PSL3; SPP; SPPL1  |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer (HP226073) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| Storage:         | Store at -20°C.   |
| Stability:       | The primer mix and template standard are stable for one year from date of shipping.   |



[View online »](#)