

## Product datasheet for **HK203754**

### hnRNP A1 (HNRNPA1) Human qPCR Template Standard (NM\_002136)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_002136
Synonyms:	ALS19; ALS20; hnRNP-A1; hnRNP A1; HNRPA1; HNRPA1L3; IBMPFD3; UP 1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230515) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)