

## Product datasheet for **HK203826**

### ERAB (HSD17B10) Human qPCR Template Standard (NM\_004493)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004493
Synonyms:	17b-HSD10; ABAD; CAMR; DUPXp11.22; ERAB; HADH2; HCD2; HSD10MD; MHBD; MRPP2; MRX17; MRX31; MRXS10; SCHAD; SDR5C1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP229833) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)