

Product datasheet for **HK203827**

HSD17B11 Human qPCR Template Standard (NM_016245)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_016245
Synonyms:	17-BETA-HSD11; 17-BETA-HSDXI; 17BHSD11; DHRS8; PAN1B; RETSDR2; SDR16C2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP212071) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)