

Product datasheet for **HK203832**

HSD17B8 Human qPCR Template Standard (NM_014234)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_014234
Synonyms:	D6S2245E; dj1033B10.9; FABG; FABGL; H2-KE6; HKE6; KE6; RING2; SDR30C1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210719) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)