

Product datasheet for **HK203839**

HSFY1 (HSFY2) Human qPCR Template Standard (NM_001001878)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001001878
Synonyms:	FLJ25453; heat shock transcription factor, Y linked 2; heat shock transcription factor 2-like; HSF2L; HSFY; HSFY, HSF2L, FLJ25453; OTTHUMP00000038975
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP232971) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)