

Product datasheet for **HK203842**

HSP70-1A (HSPA1A) Human qPCR Template Standard (NM_005345)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_005345
Synonyms:	HEL-S-103; HSP70-1; HSP70-1A; HSP70-2; HSP70.1; HSP70.2; HSP70I; HSP72; HSPA1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208499) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)