

## Product datasheet for **HK203851**

### Grp75 (HSPA9) Human qPCR Template Standard (NM\_004134)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004134
Synonyms:	CRP40; CSA; EVPLS; GRP-75; GRP75; HEL-S-124m; HSPA9B; MOT; MOT2; MTHSP75; PBP74; SAAN; SIDBA4
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207495) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)