

Product datasheet for **HK203946**

IGFBP7 Human qPCR Template Standard (NM_001553)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_001553 |
| Synonyms: | AGM; FSTL2; IBP-7; IGFBP-7; IGFBP-7v; IGFBPRP1; MAC25; PSF; RAMSVPS; TAF |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP205409) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)