

Product datasheet for **HK203973**

IL18R1 Human qPCR Template Standard (NM_003855)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003855
Synonyms:	CD218a; CDw218a; IL-1Rrp; IL-18R-alpha; IL-18Ralpha; IL1RRP; IL18RA; IL18Ralpha2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207304) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)