

Product datasheet for **HK204051**

Integrin beta 1 (ITGB1) Human qPCR Template Standard (NM_002211)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_002211
Synonyms:	CD29; FNRB; GPIIA; MDF2; MSK12; VLA-BETA; VLAB
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP232479) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)