

## Product datasheet for **HK204059**

### ITM2B Human qPCR Template Standard (NM\_021999)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_021999
Synonyms:	ABRI; BRI; BRI2; BRICD2B; E3-16; E25B; FBD; imBRI2; RDGCA
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP214303) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)