

Product datasheet for **HK204164**

KIAA0649 (PPP1R26) Human qPCR Template Standard (NM_014811)

Product data:

| | |
|-------------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_014811 |
| Synonyms: | KIAA0649; NRBE3 |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP211172) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)