

Product datasheet for **HK204326**

LGALS3BP Human qPCR Template Standard (NM_005567)

Product data:

Product Type:	qPCR Template Standards
Accession No:	NM_005567
Synonyms:	90K; BTBD17B; CyCAP; gp90; M2BP; MAC-2-BP; TANGO10B
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)