

Product datasheet for **HK204522**

MAGEB2 Human qPCR Template Standard (NM_002364)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_002364
Synonyms:	CT3.2; DAM6; MAGE-XP-2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206065) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)