

## Product datasheet for **HK204555**

### MEK4 (MAP2K4) Human qPCR Template Standard (NM\_003010)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003010
Synonyms:	JNKK; JNKK1; MAPKK4; MEK4; MKK4; PRKMK4; SAPKK-1; SAPKK1; SEK1; SERK1; SKK1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206611) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)