

Product datasheet for **HK204594**

Methionyl tRNA synthetase (MARS) Human qPCR Template Standard (NM_004990)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004990
Synonyms:	CMT2U; ILFS2; ILLD; MARS; METRS; MRS; MTRNS; SPG70
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208215) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)