

Product datasheet for **HK204631**

Mucolipin 1 (MCOLN1) Human qPCR Template Standard (NM_020533)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_020533
Synonyms:	MG-2; ML1; ML4; MLIV; MST080; MSTP080; TRP-ML1; TRPM-L1; TRPML1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213716) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)