

Product datasheet for **HK204784**

MMP15 Human qPCR Template Standard (NM_002428)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_002428
Synonyms:	MMP-15; MT2-MMP; MT2MMP; MTMMP2; SMCP-2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206114) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)