

Product datasheet for **HK204894**

MRPS26 Human qPCR Template Standard (NM_030811)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_030811
Synonyms:	C20orf193; dj534B8.3; GI008; MRP-S13; MRP-S26; MRPS13; NY-BR-87; RPMS13
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP215431) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)