

Product datasheet for **HK204897**

MRPS28 Human qPCR Template Standard (NM_014018)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_014018
Synonyms:	COXPD47; HSPC007; MRP-S28; MRP-S35; MRPS35
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210586) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)