

Product datasheet for **HK205117**

NFKB1 Human qPCR Template Standard (NM_003998)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003998
Synonyms:	CVID12; EBP-1; KBF1; NF-kappa-B1; NF-kappaB; NF-kappabeta; NF-kB; NF-kB1; NFkappaB; NFKB-p50; NFKB-p105
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207409) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)