

Product datasheet for **HK205192**

NPY5R Human qPCR Template Standard (NM_006174)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_006174
Synonyms:	NPY5-R; NPYR5; NPYY5-R
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209161) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)