

Product datasheet for **HK205209**

NUR77 (NR4A1) Human qPCR Template Standard (NM_173157)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_173157
Synonyms:	GFRP1; HMR; N10; NAK-1; NGFIB; NP10; NUR77; TR3
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205874) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)