

Product datasheet for **HK205309**

Sigma1 receptor (SIGMAR1) Human qPCR Template Standard (NM_005866)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_005866
Synonyms:	ALS16; DSMA2; hSigmaR1; OPRS1; SIG-1R; sigma1R; SR-BP; SR-BP1; SRBP
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208924) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)