

## Product datasheet for **HK205362**

### **P2Y10 (P2RY10) Human qPCR Template Standard (NM\_014499)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_014499
<b>Synonyms:</b>	LYPSR2; P2Y10
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP210935) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)