

## Product datasheet for **HK205397**

### PAK4 Human qPCR Template Standard (NM\_001014833)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001014833
Synonyms:	p21 protein (Cdc42/Rac)-activated kinase 4; p21(CDKN1A)-activated kinase 4; p21-activated kinase 4; protein kinase related to <i>S. cerevisiae</i> STE20, effector for Cdc42Hs
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP227807) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)