

## Product datasheet for **HK205473**

### PD L2 (PDCD1LG2) Human qPCR Template Standard (NM\_025239)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_025239
Synonyms:	B7DC; bA574F11.2; Btdc; CD273; PD-L2; PDCD1L2; PDL2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP215316) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)