

## Product datasheet for **HK205583**

### PHKA2 Human qPCR Template Standard (NM\_000292)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000292
Synonyms:	GSD9A; PHK; PYK; PYKL; XLG; XLG2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200273) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)