

## Product datasheet for **HK205653**

### Calcium independent Phospholipase A2 (PLA2G6) Human qPCR Template Standard (NM\_003560)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003560
Synonyms:	CaI-PLA2; GVI; INAD1; iPLA2; IPLA2-VIA; iPLA2beta; NBIA2; NBIA2A; NBIA2B; PARK14; PLA2; PNPLA9
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230089) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)