

## Product datasheet for **HK205704**

### Probable hydrolase PNKD (PNKD) Human qPCR Template Standard (NM\_015488)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_015488
Synonyms:	BRP17; DYT8; FKSG19; FPD1; KIPP1184; MR-1; MR-1S; MR1; PDC; PKND1; PNKD1; R1; TAHCCP2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230623) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)