

## Product datasheet for **HK205885**

### PRMT5 Human qPCR Template Standard (NM\_006109)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_006109
Synonyms:	HRMT1L5; HSL7; IBP72; JBP1; SKB1; SKB1Hs
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230776) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)