

Product datasheet for **HK206024**

SHP1 (PTPN6) Human qPCR Template Standard (NM_080548)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_080548
Synonyms:	HCP; HCPH; HPTP1C; PTP-1C; SH-PTP1; SHP-1; SHP-1L; SHP1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP231470) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)