

Product datasheet for **HK206153**

Epac2 (RAPGEF4) Human qPCR Template Standard (NM_007023)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_007023
Synonyms:	CAMP-GEFII; CGEF2; EPAC; EPAC 2; EPAC2; Nbla00496
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230143) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)