

Product datasheet for **HK206170**

Rb (RB1) Human qPCR Template Standard (NM_000321)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000321
Synonyms:	OSRC; p105-Rb; p110-RB1; pp110; PPP1R130; pRb; RB
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200302) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)