

Product datasheet for **HK206214**

RDH11 Human qPCR Template Standard (NM_016026)

Product data:

| | |
|-------------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_016026 |
| Synonyms: | ARSDR1; CGI82; HCBP12; MDT1; PSDR1; RALR1; RDJCSS; SCALD; SDR7C1 |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP211921) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)