

Product datasheet for **HK206486**

TIP49A (RUVBL1) Human qPCR Template Standard (NM_003707)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003707
Synonyms:	ECP-54; ECP54; INO80H; NMP 238; NMP238; PONTIN; Pontin52; RVBL1; TIH1; TIP49; TIP49A
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207178) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

