

## Product datasheet for **HK206540**

### SCARB1 Human qPCR Template Standard (NM\_005505)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_005505
Synonyms:	CD36L1; CLA-1; CLA1; HDLQTL6; SR-BI; SRB1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230174) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)