

Product datasheet for **HK206719**

Melanoma gp100 (PMEL) Human qPCR Template Standard (NM_006928)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_006928
Synonyms:	D12S53E; gp100; ME20; ME20-M; ME20M; P1; P100; PMEL17; SI; SIL; SILV
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209772) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)