

Product datasheet for **HK206762**

Solute carrier family 22 member 18 (SLC22A18) Human qPCR Template Standard (NM_183233)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_183233
Synonyms:	BWR1A; BWSCR1A; HET; IMPT1; ITM; ORCTL2; p45-BWR1A; SLC22A1L; TSSC5
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206215) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)