

Product datasheet for **HK206785**

Adenine Nucleotide Translocator 2 (SLC25A5) Human qPCR Template Standard (NM_001152)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001152
Synonyms:	2F1; AAC2; ANT2; T2; T3
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205090) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)