

## Product datasheet for **HK206883**

### SMCR7L (MIEF1) Human qPCR Template Standard (NM\_019008)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_019008
Synonyms:	AltMIEF1; dJ1104E15.3; HSU79252; MID51; MIEF1-MP; SMCR7L
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213340) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)