

Product datasheet for **HK206910**

SNRP70 (SNRNP70) Human qPCR Template Standard (NM_003089)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003089
Synonyms:	RNPUIZ; RPU1; Snp1; SNRP70; U1-70K; U1AP; U1RNP; U170K
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP233850) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

