

Product datasheet for **HK207068**

CDw75 (ST6GAL1) Human qPCR Template Standard (NM_003032)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_003032 |
| Synonyms: | SIAT1; ST6Gall; ST6N |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP232459) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)