

Product datasheet for **HK207278**

TEX11 Human qPCR Template Standard (NM_031276)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_031276
Synonyms:	MZIP4; SPGFX2; Spo22; TGC1; TSGA3; ZIP4; ZIP4H
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP201265) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)