

## Product datasheet for **HK207348**

### Transketolase (TKT) Human qPCR Template Standard (NM\_001064)

#### Product data:

|                         |                                                                                                                                                                                                                                                                                               |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>    | qPCR Template Standards                                                                                                                                                                                                                                                                       |
| <b>ACCN:</b>            | NM_001064                                                                                                                                                                                                                                                                                     |
| <b>Synonyms:</b>        | HEL-S-48; HEL107; SDDHD; TK; TKT1                                                                                                                                                                                                                                                             |
| <b>Components:</b>      | 1 vial of lyophilized SybGREEN qPCR primer (HP203785) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| <b>Quality Control:</b> | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT                                                                                                                                                                                                                      |
| <b>Storage:</b>         | Store at -20°C.                                                                                                                                                                                                                                                                               |
| <b>Stability:</b>       | The primer mix and template standard are stable for one year from date of shipping.                                                                                                                                                                                                           |

