

Product datasheet for **HK207527**

PRPK (TP53RK) Human qPCR Template Standard (NM_033550)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_033550
Synonyms:	BUD32; C20orf64; dj101A2; GAMOS4; Nori-2; Nori-2p; PRPK; TPRKB
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP216465) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)