

Product datasheet for **HK207625**

TM4SF2 (TSPAN7) Human qPCR Template Standard (NM_004615)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004615
Synonyms:	A15; CCG-B7; CD231; DXS1692E; MRX58; MXS1; TALLA-1; TM4SF2; TM4SF2b
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207903) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)