

Product datasheet for **HK207654**

ETS1 associated protein II (TDP2) Human qPCR Template Standard (NM_016614)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_016614
Synonyms:	AD022; dj30M3.3; EAP2; EAPII; hTDP2; TTRAP
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP212301) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)