

Product datasheet for **HK207661**

beta V Tubulin (TUBB) Human qPCR Template Standard (NM_178014)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_178014
Synonyms:	CDCBM6; CSCSC1; M40; OK/SW-cl.56; TUBB1; TUBB5
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP218747) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)